


Seed Respiration Analyser (SRA)

Germination and Respiration Data Analysis

The SRA can support traditional seed testing, providing information unavailable through traditional means. The SRA has four discrete temperature zones, making it possible to easily compare germination and seed respiration responses to multiple temperature regimes. Easy to use SRA data analysis tools are freely available from UC Davis, Department of Plant sciences. The analysis software tool delivers automatically a results report as shown below.




O₂ depletion curves. Each line represents one single seed. The dashed lines are the average (black) and median (red) curves. After an initial period of low O₂ consumption, the seeds accelerate respiration (steep curve). The heterogeneity of the batch can be clearly seen.

Experiment Information

Description Run: **02-Mar-21**
 Description Plate: **Cress**
 ID Run: **52** Plate #: **1** Plate type: **35**
 Imbibition date: **3-2-21 9:30 AM**
 Starting date: **3-2-21 11:30 AM**
 Finishing date: **3-6-21 5:30 PM**
 Temperature: **25** Lot #: **0**


Variety	Lot 1		
Treatment 1	Description: Agar	Value: 400	
Treatment 2	Description: Temp	Value: 25	
Additional Info			



	x min	x max
SMR	5	15
OMR	45	60


SMR	0.24	%/hr
IMT	36.91	hr
OMR	1.62	%/hr
R50	62.21	hr
AUC50	53.36	%.hr
RGT	92.89	hr

Average O₂ depletion curve (dashed line) with Standard Error of the Mean. In addition, several parameters derived from the O₂ depletion curve are given in the table (see glossary for explanation). SMR and OMR (start and maximal metabolism rate) are shown in the graph as well.



Results report

R75 POD Curve - Values					
POD R75(50)	46.50	STDEV		T90-10	35.25
Rate POD R75(50)	2.15%			Final (%)	97%
R75	52.12	15.85		AUC75	45.32
R50 POD Curve - Values					
POD R50(50)	59.50	STDEV		T90-10	44.50
Rate R50(50)	1.68%			Final (%)	91%
R50	62.16	11.892		AUC50	51.58
R25 POD Curve - Values					
POD R25(50)	75.50	STDEV		T90-10	#N/A
Rate POD R25(50)	1.32%			Final (%)	89%
R25	74.97	10.94		AUC25	55.79



Population O₂ Depletion (POD) curves showing the accumulation of the % of seeds reaching 75, 50 and 25 % O₂ level in the vial (from the O₂ single seed depletion curves). Step lines (close together) represent high homogeneity. Table on the right shows different parameter values and statistical information from the POD analysis

Glossary

- SMR - Starting Metabolism Rate
- IMT - Increased Metabolism Time
- OMR - Oxygen Metabolism Rate
- RGT - Relative Germination Time
- R75(50)/R50(50)/R25(50) - Time to 50% of seeds reaching 75/50/25% oxygen level
- T90-10 - Time between 10% and 90% of seeds reaching 75/50/25% oxygen level
- Rate to R75(50)/R50(50)/R25(50) - Inverse of time to 50% of seeds reaching 75/50/25% oxygen level
- R75/50/25 - Individual-seed time to reach 75/50/25% oxygen level
- AUC75/50/25 - Individual-seed area under the curve to 75/50/25% oxygen level